



# Filtering Purification

## Aire~LifeLine™

**In-line filter for cylinder and house gasses  
(O<sub>2</sub>, N<sub>2</sub>, CO<sub>2</sub>, and tri-gas)**

- LifeAire's advanced purification technology inside
- Minimum 6 month filtering capacity
- Multi-step filtration for added protection
- Delivers broadest VOC and particulate protection
- Replacement only twice a year
- Quick connects into existing set-ups
- Off-gassed in advance
- Magnetic and wall-mounting bracket option
- CE certification

For more information, contact [emarsh@lifeaire.com](mailto:emarsh@lifeaire.com)  
Visit [www.lifeaire.com](http://www.lifeaire.com) or call 484-224-3042

LifeAire Systems